


USER MANUAL

RED LIGHT THERAPY PANEL



CONTENTS

I	Product Parameters	01
II	Control Panel Operating Instructions	02
III	Accessory Instruction	06
IV	Important Usage Guidelines	07



○ Product Parameters

Item	T70-COB
LED Power	360 W
LED QTY	68+COB
Input Voltage	100-240AC
Frequency	50/60HZ
Lens	30°
Irradiance	0inch ≥ 200mW/cm ² ±10%
Dimension	352*235*70 mm
Spectrum	485nm+630nm+660nm+727nm +810nm+850nm+940nm
Actual Power	110 W±10%


○ Control Panel Operating Instructions

Power On Page



This page will be displayed when you start up, and after a short time you will go to the home page

Home Page



Display remaining time and enter the timed page directly by single-clicking

Display the current value

Start / Pause

Click to enter the *Set* page

Click to enter the *Mode* page

Mode Page



Fitness mode



Good build mode



Muscle recovery mode



Inflammation mode



Skin care mode



Sleep mode



Growth hair mode



Pet mode



Turn on/off the connection

M₁

Custom mode 1

M₂

Custom mode 2



Return

Click the icon to enter the corresponding mode and automatically jump to the **home page**

Set Page



Click to enter the **Timer page**

Click to enter the **Pulse page**

Click to enter the **Dimming page**

Save the current value (brightness, frequency) as a custom mode, a total of two modes can be saved, saving the new mode will replace the old mode

Click the icon to jump to the corresponding page, Press and hold "⚙️" for 6 seconds to turn on/off the sound

⚙️ **Set page**


🕒 **Timer page**

🏠 **Home page**

📶 **Pulse page**

☀️ **Dimming page**

Timer Page



Displays the current remaining time

Drag the white slider to adjust the time, up to 30 minutes. Click the + - button on both sides to adjust the time in units of 1 minute

Click the button for shortcut Settings

Click the icon to jump to the corresponding page,
Press and hold " ⚙ " for 6 seconds to turn on/off the sound

⚙ *Set page* ⌚ *Timer page* 🏠 *Home page* 📶 *Pulse page* ☀ *Dimming Page*

Dimming Page



Drag the slider to adjust the brightness

Drag the slider to adjust the brightness

Drag the slider to adjust the brightness

Drag the slider to adjust the brightness

Click the button for a quick setup (What changed was the overall brightness.)

Click the icon to jump to the corresponding page,
Press and hold " ⚙ " for 6 seconds to turn on/off the sound

⚙ *Set page* ⌚ *Timer page* 🏠 *Home page* 📶 *Pulse page* ☀ *Dimming Page*

Pulse Page




Click the button to turn on/off the pulse, the current pulse frequency will be displayed in the box, click the box can jump to the Number page, set the pulse frequency freely in the digital page

Click the button to quickly set the parameters (what changes is the overall brightness.)

Click the icon to jump to the corresponding page,
Press and hold "⚙️" for 6 seconds to turn on/off the sound

⚙️ *Set page* ⌚ *Timer page* 🏠 *Home page* 📶 *Pulse page* ☀️ *Dimming Page*

Digital Page



Click the number button to set the frequency, the value will be displayed in the box, click OK button, the setting takes effect, and the pulse page will be automatically returned. Click the arrow button to go back a number. Click the House button to unset the return pulse page.

The maximum value is 300. Input number exceeding 300 will automatically change to 300.

Accessory Instruction

Hanging rope 30cm



Step 1


Align the smaller end of the double-ended screw with the reserved screw hole on the red light treatment panel and screw into it.

Step 2

Screw the sling into the secured double-headed screws.

Step 3

It can be hung anywhere.



Door hook




Step 1

Attach the hook to the door.

Step 2

Hang the red light physiotherapy instrument with the hanging rope on the door hook.



Halyard 200cm



Step 1


Connect the lower silver hook to the red light treatment panel sling.

Step 2

Pull the spring arm on the side of the device upward to adjust the length of the rope to the desired size.

Step 3

The top silver hook can be hung anywhere.



○ Important Usage Guidelines

Power Connection:

Ensure the power cable is securely connected to both the light and the power socket. Switch on the light using the power button. Use the buttons on the control panel or the remote controller to operate the light. Before use, carefully read the manual to understand the functions of the control panel.

● **Distance and Time Recommendations.**

Optimal results are achieved with sessions at a distance of 15 to 30 cm and for 10 to 30 minutes per treated area.

Suggested session duration: 10 to 20 minutes per side (e.g., front and back).

Recommended frequency: once a day, if possible.

Adjust time and distance based on individual responses to the therapy.

● **Eye Protection**

We recommend wearing the included darkened glasses provided in the set for optimal eyeprotection.

● **Clothing**

Ensure the treated area is fully exposed to the light and free from clothing or fabric.

You may cover non-treated areas for comfort though it is not necessary.

● **Safety**

Red Light Therapy is safe, non-invasive, and free from chemicals, medicines, and UV radiation.

It is not associated with harmful side effects, making it a safe alternative to conventional treatments.

Following these guidelines will help you achieve the best and safest results.

Troubleshooting

Below are solutions to common issues you may encounter with the panel:

1. Why is my panel not lighting up?

Ensure the power cord is fully inserted into the connector. Check if the other end of the power cord is securely connected to the wall socket. Test the light by plugging it into a different wall socket. Note: The red light is visible, but the infrared light is invisible to the human eye.

2. Why do some of the LEDs look like they are not working?

All lights use NIR (Near Infrared) LEDs, which may appear dim or non-functional because infrared light is virtually invisible to the human eye.

Although they may not look active, NR LEDs are incredibly bright and effective for your body. This appearance is completely normal and does not indicate any malfunction.

Safety Information

● Assembly

Ensure the lamp is assembled with the help of at least two people to prevent falls or damage. One person should lift the lamp while the other aligns and tightens the bolts securely.

● Safety Considerations

described in this manual will void the warranty and could lead to serious injury. As the panel is an electrical device, adhere to the following safety precautions to avoid electric shock:

- Do not expose the device to water or submerge it in any liquid.
- Avoid washing electrical components with water or other liquids.
- Never touch the device, plugs, or switches with wet hands.
- Do not apply or operate the device on wet or damp skin.
- Avoid pulling, carrying, or lifting the device by its power cord. If the cord is damaged, discontinue use immediately.
- Do not use the device if it has been dropped, left outdoors, or exposed to water or liquids.
- Keep the device away from heated surfaces.
- Do not modify or attempt to repair the device. This equipment contains no user-serviceable parts.
- Avoid using or storing the device in environments where temperatures exceed 130°F(54°C)

THANKS